

Human EGFR (erbb1) protein (Recombinant) (STJP000302)

STJP000302

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human EGFR (erbb1)-protein was developed from insect cells. For use in research applications.
Host/Source	Insect cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 M unfiltered 20mM Tris, 50mM NaCl, pH 7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. It is recommended to add 0.1% HSA or BSA for long term storage. Avoid repeated freeze/thaw cycles. NA
Instruction	

TARGET INFORMATION

Gene ID	1956
Gene Symbol	EGFR
Uniprot ID	EGFR_HUMAN
Immunogen Sequence	