

Human Platelet Derived Growth Factor-BB protein (Recombinant) (STJP000308)

STJP000308

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Platelet Derived Growth Factor-BB-protein was developed from s.cerevisiae. For use in research applications.
Host/Source	S.cerevisiae

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered.02M acetic acid solution pH.5.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC. NA
Storage	Can be stored in working aliquots at°C-°C C for one month, or at-20°C C for six months, with a carrier protein without detectable
Instruction	loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	5155
Gene Symbol	PDGFB
Uniprot ID	PDGFB_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081