

Human TGF Beta-1 protein (Recombinant) (Fc-Tag) (STJP000315)

STJP000315

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human TGF Beta-1-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered 100mM glycine PH 3.0.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC. NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months. Avoid repeated freeze/thaw cycles.
Instruction	NA

TARGET INFORMATION

Gene ID	7040
Gene Symbol	TGFB1
Uniprot ID	TGFB1_HUMAN
Immunogen	
Sequence	