

Mouse TGFB1 protein (Recombinant) (STJP000350)

STJP000350

GENERAL INFORMATION



Product Type	Proteins
Short Description	Recombinant-Mouse TGFB1-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Gene ID	21803
Gene Symbol	Tgfb1
Uniprot ID	TGFB1_MOUSE
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081