

Human CD10 (Neprilysin) protein (Recombinant) (STJP000386)

STJP000386

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD10 (Neprilysin)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	0.2 Mu m filtered solution in 20mM Tris, 100mM NaCL, pH8.0 with 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months, with 0.1% HSA as carrier protein without
Instruction	detectable loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	4311
Gene Symbol	MME
Uniprot ID	NEP_HUMAN
Immunogen Sequence	