

Human CD32a (FcγRIIa) protein (Recombinant) (STJP000413)

STJP000413

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD32a (FcγRIIa)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised 0.2 µg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	Lyophilized CD32a, as supplied, should be stored at -20°C to -70°C. Upon reconstitution, working aliquots may be stored at 2°C-8°C.
Instruction	for one month, or six months at -20°C without detectable loss of activity. It is recommended to add 0.1% BSA for long term.

TARGET INFORMATION

Gene ID	2212
Gene Symbol	FCGR2A
Uniprot ID	FCG2A_HUMAN
Immunogen	
Sequence	