

## Human Hepatocyte Growth Factor protein (Recombinant) (STJP000426)

STJP000426

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human Hepatocyte Growth Factor-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.2 Mu m filtered solution in 0.5% glycine, 0.5% sucrose,.01% Tween80, mM Glutamic acid, pH7.5.
<b>Purification</b>	
<b>Dilution Range</b>	>97%, as determined by SDS-PAGE and HPLC NA
<b>Storage</b>	Can be stored in working aliquots at°C-°C C for one month, or at-20°C C for six months, with a carrier protein without detectable
<b>Instruction</b>	loss of activity. Avoid repeated freeze/thaw cycles. NA

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3082</a>
<b>Gene Symbol</b>	<a href="#">HGF</a>
<b>Uniprot ID</b>	<a href="#">HGF_HUMAN</a>
<b>Immunogen</b>	
<b>Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081