

Human Follistatin protein (Recombinant) (STJP000469)

STJP000469

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Follistatin-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised at 1mg/ml containing no additives.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months, with a carrier protein without detectable
Instruction	loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	10468
Gene Symbol	FST
Uniprot ID	FST_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081