

Mouse Stem Cell Factor protein (Recombinant) (STJP000479) STJP000479

GENERAL INFORMATION

Product Type Proteins

 Short
 Recombinant-Mouse Stem Cell Factor-protein was developed from insect cells. For use in research applications.

 Description
 Insect cells

PRODUCT PROPERTIES

Concentration Formulation

Formulation Lyophilised from 0.2 Mu m filtered solution in mM Tris, % trehalose, pH7.5. Purification

 Purification

 Dilution Range

 Storage

 Can be stored in working aliquots at°C-°C C for one month, or at-20°C C for six months, with a carrier protein without detectable

Instruction loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes. St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081