

Rat Lipocalin-2 protein (Recombinant) (STJP000508)

STJP000508

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | Recombinant-Rat Lipocalin-2-protein was developed from hek293. For use in research applications. |
| Host/Source | HEK293 |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | Lyophilised from 0.2 µm filtered PBS solution, pH7.2, 5% Trehalose. |
| Purification | |
| Dilution Range | >95%, as determined by SDS-PAGE and HPLC NA |
| Storage Instruction | Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA |

TARGET INFORMATION

| | |
|--------------------|----------|
| Gene ID | 170496 |
| Gene Symbol | Lcn2 |
| Uniprot ID | NGAL_RAT |
| Immunogen | |
| Sequence | |

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081