

## Mouse VEGFR2 protein (Recombinant) (Fc-Tag) (STJP000524)

STJP000524

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse VEGFR2-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC NA
<b>Storage Instruction</b>	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">16542</a>
<b>Gene Symbol</b>	<a href="#">Kdr</a>
<b>Uniprot ID</b>	<a href="#">VGFR2_MOUSE</a>
<b>Immunogen</b>	
<b>Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081