

Human TGF Beta-1 receptor 2 (TGRBR2) protein (Recombinant) (STJP000530)

STJP000530

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human TGF Beta-1 receptor 2 (TGRBR2)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 M filtered 100mM Glycine, 10mM NaCl, 50mM Tris, pH 7.5.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC. NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C for one month, or at -20°C for six months. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	7048
Gene Symbol	TGFB2
Uniprot ID	TGFR2_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081