

**Mouse CD27 (TNFRSF7) protein (Recombinant) (STJP000538)**

STJP000538

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | Recombinant-Mouse CD27 (TNFRSF7)-protein was developed from hek293. For use in research applications. |
| <b>Host/Source</b>       | HEK293  |

**PRODUCT PROPERTIES**

|                            |  |
|----------------------------|--|
| <b>Concentration</b>       |  |
| <b>Formulation</b>         | Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      | >97%, as determined by SDS-PAGE and HPLC NA  |
| <b>Storage Instruction</b> | Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA |

**TARGET INFORMATION**

|                    |                            |
|--------------------|----------------------------|
| <b>Gene ID</b>     | <a href="#">21940</a>      |
| <b>Gene Symbol</b> | <a href="#">Cd27</a>       |
| <b>Uniprot ID</b>  | <a href="#">CD27_MOUSE</a> |
| <b>Immunogen</b>   |                            |
| <b>Sequence</b>    |                            |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081