

SARS Spike protein S1 protein (Recombinant) (STJP000565)

STJP000565

GENERAL INFORMATION


Product Type	Proteins
Short Description	Recombinant-SARS Spike protein S1-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES


Concentration	
Formulation	
Purification	
Dilution Range	>95%, as determined by SDS-PAGE. NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence



Tris-Bis PAGE under reducing condition. The purity is greater than 95%



The purity of SARS hCov Spike S1 is greater than 95% as determined by SEC-HPLC.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081