

Rat IL-1 receptor 1 (CD121a) protein (Recombinant) (STJP000591)

STJP000591

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Rat IL-1 receptor 1 (CD121a)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC NA
Storage Instruction	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	
Gene Symbol	Il1r1
Uniprot ID	IL1R1_RAT
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081