

## Hemagglutinin Influenza A-Indonesia-05-2005-H5N1 protein (Recombinant) (STJP000616)

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Hemagglutinin Influenza A-Indonesia-05-2005-H5N1-protein was developed from insect cells. For use in research applications.
<b>Host/Source</b>	Insect cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from PBS pH7.4, 5% trehalose, 5% mannitol.
<b>Purification</b>	
<b>Dilution Range</b>	>90%, as determined by SDS-PAGE and HPLC NA
<b>Storage</b>	Can be stored in working aliquots at 2°C-8°C C for one month, or at-20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Sequence</b>