

Hemagglutinin Influenza A-Vietnam 1203-2004-H5N1 protein (Recombinant) (STJP000622)

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Hemagglutinin Influenza A-Vietnam 1203-2004-H5N1 -protein was developed from insect cells. For use in research applications.
Host/Source	Insect cells

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from PBS pH7.4, 5% trehalose, 5% mannitol.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles. NA
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence