

Mouse IL-10 receptor beta (IL10RB) protein (Recombinant) (STJP000628)

STJP000628

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-10 receptor beta (IL10RB)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 Mu m filtered PBS solution pH7.4.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC NA
Storage	Can be stored in working aliquots at°C-°C C for one month, or at-20°C C for six months, with a carrier protein without detectable
Instruction	loss of activity. Avoid repeated freeze/thaw cycles. NA

TARGET INFORMATION

Gene ID	
Gene Symbol	Il10rb
Uniprot ID	I10R2_MOUSE
Immunogen	
Sequence	