

**Human CD114 (GCSFR) protein (Recombinant) (STJP000675)**

STJP000675

**GENERAL INFORMATION**

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human CD114 (GCSFR)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

**PRODUCT PROPERTIES**

<b>Concentration</b>	
<b>Formulation</b>	0.2 Mu m filtered PBS solution, pH7.2.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC NA
<b>Storage Instruction</b>	Can be stored in working aliquots at 2°C-8°C C for one month, or at -20°C to -70°C for 1 year. Avoid repeated freeze/thaw cycles. NA

**TARGET INFORMATION**

<b>Gene ID</b>	1441
<b>Gene Symbol</b>	CSF3R
<b>Uniprot ID</b>	CSF3R_HUMAN
<b>Immunogen</b>	
<b>Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081