

Human EPCR protein (Recombinant) (Fc-Tag) (STJP000707)

STJP000707

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human EPCR-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC NA
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C. NA

TARGET INFORMATION

Gene ID	10544
Gene Symbol	PROCR
Uniprot ID	EPCR_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081