

**Human VEGF-A/VEGF121 (K321N) protein (Recombinant) (No tag) (STJP000791)**  
STJP000791

**GENERAL INFORMATION**





**Product Type** Proteins  
**Short Description**  
**Host/Source** E. coli

**PRODUCT PROPERTIES**

**Concentration**  
**Formulation** Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.  
**Purification**  
**Dilution Range** > 97% by SDS-PAGE.  
**Storage Instruction** Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

**Gene ID** 7422  
**Gene Symbol** VEGFA  
**Uniprot ID** VEGFA\_HUMAN  
**Immunogen Sequence**



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081