

Human IDO-1 protein (Recombinant) (No tag) (STJP000797)

STJP000797

GENERAL INFORMATION

Product Type Proteins
Short Description
Host/Source E. coli

PRODUCT PROPERTIES

Concentration

Formulation Lyophilised from a 0.22 Mu m filtered solution of 50mM Tris, 150mM NaCl, pH 8.0.


Purification

Dilution Range > 97% by SDS-PAGE.

Storage Instruction Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 3620
Gene Symbol IDO1
Uniprot ID 123O1_HUMAN
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes. St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081