

## Human ERK2/MAPK1 protein (Recombinant) (N-GST) (STJP000800)

STJP000800

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	
<b>Host/Source</b>	E. coli

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from a 0.22 µg/ml filtered solution of 50mM Tris, 150mM NaCl, 0.1mM EDTA, 0.25mM DTT, 0.1mM AEBSF, pH 8.0.
<b>Purification</b>	
<b>Dilution Range</b>	> 70% by SDS-PAGE.
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	5594
<b>Gene Symbol</b>	MAPK1
<b>Uniprot ID</b>	MK01_HUMAN
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081