

Human IL-11 protein (Active Recombinant) (No tag) (STJP000825)
STJP000825

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source E. coli

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Purification
Dilution Range > 97% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 3589
Gene Symbol IL11
Uniprot ID IL11_HUMAN
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081