

**Human FKBP1A/FKBP12 protein (Recombinant) (No tag) (STJP000826)**

STJP000826

**GENERAL INFORMATION**


|                          |          |
|--------------------------|----------|
| <b>Product Type</b>      | Proteins |
| <b>Short Description</b> |          |
| <b>Host/Source</b>       | E. coli  |

**PRODUCT PROPERTIES**

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Concentration</b>       |                                                                                                |
| <b>Formulation</b>         | Lyophilised from a 0.22 M filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.                  |
| <b>Purification</b>        |                                                                                                |
| <b>Dilution Range</b>      | > 90% by SDS-PAGE.                                                                             |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |             |
|---------------------------|-------------|
| <b>Gene ID</b>            | 2280        |
| <b>Gene Symbol</b>        | FKBP1A      |
| <b>Uniprot ID</b>         | FKB1A_HUMAN |
| <b>Immunogen Sequence</b> |             |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081