

Human Alpha-2-MRAP/LRPAP1 protein (Recombinant) (C-His) (STJP000833)

STJP000833

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E. coli

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from a 0.22 M filtered solution of 20mM Tris, 150mM NaCl, 5% Glycerol, pH 8.0.
Purification	
Dilution Range	> 87% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4043
Gene Symbol	LRPAP1
Uniprot ID	AMRP_HUMAN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081