

Human GDNFR-alpha-1/GFRA1 protein (Recombinant) (C-hFc&His) (STJP000878)
STJP000878

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 80% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 2674
Gene Symbol GFRA1
Uniprot ID GFRA1_HUMAN
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081