

**Human uPA/PLAU protein (Recombinant) (C-His) (STJP000921)**  
STJP000921

**GENERAL INFORMATION**





**Product Type** Proteins  
**Short Description**  
**Host/Source** HEK293 cells

**PRODUCT PROPERTIES**

**Concentration**  
**Formulation** Supplied as a 0.22  $\mu$ g/ml filtered solution in 20mM HEPES, 150mM NaCl, 2mM CaCl<sub>2</sub>, 10% Glycerol, pH 7.5  
**Purification**  
**Dilution Range** > 97% by SDS-PAGE.  
**Storage Instruction** Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

**Gene ID** 5328  
**Gene Symbol** PLAU  
**Uniprot ID** UROK\_HUMAN  
**Immunogen Sequence**



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081