

Human Cathepsin S protein (Active Recombinant) (C-His) (STJP000955)

STJP000955

GENERAL INFORMATION


| | |
|--------------------------|--------------|
| Product Type | Proteins |
| Short Description | |
| Host/Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Concentration | |
| Formulation | Lyophilised from a 0.22 µg/ml filtered solution of 50mM MES, 100mM NaCl, pH 6.5. |
| Purification | |
| Dilution Range | > 95% by SDS-PAGE. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 1520 |
| Gene Symbol | CTSS |
| Uniprot ID | CATS_HUMAN |
| Immunogen Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081