

## Human TNFRSF3/TNFR-III/LTBR protein (Recombinant) (C-hFc&His) (STJP001190)

STJP001190

### GENERAL INFORMATION





|                          |              |
|--------------------------|--------------|
| <b>Product Type</b>      | Proteins     |
| <b>Short Description</b> |              |
| <b>Host/Source</b>       | HEK293 cells |

### PRODUCT PROPERTIES

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| <b>Concentration</b>       |                                                                                                |
| <b>Formulation</b>         | Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.                                    |
| <b>Purification</b>        |                                                                                                |
| <b>Dilution Range</b>      | > 95% by SDS-PAGE.                                                                             |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            | 4055       |
| <b>Gene Symbol</b>        | LTBR       |
| <b>Uniprot ID</b>         | TNR3_HUMAN |
| <b>Immunogen Sequence</b> |            |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081