

Human GTPase KRas/KRAS (G12V) protein (Recombinant) (N-His) (STJP001252)

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source E. coli

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Purification
Dilution Range > 85% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 3845
Gene Symbol KRAS
Uniprot ID RASK_HUMAN
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081