

Human TET2 protein (Recombinant) (Please contact us for more information.) (STJP001260)

STJP001260

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	HEK293 cells

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.22 μ g/ml filtered solution in 50mM HEPES, 300mM NaCl, 1mM TCEP, 10% Glycerol, pH 8.0.
Purification	
Dilution Range	> 90% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54790
Gene Symbol	TET2
Uniprot ID	TET2_HUMAN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081