

SARS-CoV-2 3C-likease protein (Recombinant) (N-His&Avi) (STJP001273)

STJP001273

GENERAL INFORMATION


| | |
|--------------------------|----------|
| Product Type | Proteins |
| Short Description | |
| Host/Source | E. coli |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Concentration | |
| Formulation | Supplied as a 0.22 μ g/ml filtered solution in 20mM HEPES, 120mM NaCl, 10% Glycerol, 2mM TCEP, pH 7.4 |
| Purification | |
| Dilution Range | > 95% by SDS-PAGE. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081