

**SARS-CoV-2 Spike S1+S2 ECD (S-ECD) (D614G) protein (Active Recombinant) (C-His)**  
(STJP001287)  
STJP001287

**GENERAL INFORMATION**





**Product Type** Proteins  
**Short Description**  
**Host/Source** HEK293 cells

**PRODUCT PROPERTIES**

**Concentration**  
**Formulation** Supplied as a 0.22  $\mu$ M filtered solution in PBS, pH 7.4.  
**Purification**  
**Dilution Range** > 95% by SDS-PAGE.  
**Storage Instruction** Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen Sequence**



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081