

MERS-CoV Nucleocapsid protein (Recombinant) (C-His) (STJP001291)

STJP001291

GENERAL INFORMATION


Product Type	Proteins
Short Description	
Host/Source	E. coli

PRODUCT PROPERTIES

Concentration	
Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris, 500mM NaCl, 0.1mM EDTA, 20% Glycerol, pH 8.0.
Purification	
Dilution Range	> 90% by SDS-PAGE.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081