

Mouse Dkk-1 protein (Recombinant) (C-His) (STJP001341)
STJP001341

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Supplied as a 0.22 μ M filtered solution in 20mM Tris, 500mM NaCl, 10% Glycerol, pH 8.0.
Purification
Dilution Range > 90% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 13380
Gene Symbol Dkk1
Uniprot ID DKK1_MOUSE
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081