

Human Mature TGF-beta 1 protein (Active Recombinant) (No tag) (STJP001449)
STJP001449

GENERAL INFORMATION


Product Type Proteins
Short Description
Host/Source HEK293 cells

PRODUCT PROPERTIES

Concentration
Formulation Lyophilised from a 0.22 Mu m filtered solution of 50 mM Glycine, 150 mM NaCl, pH 3.5.
Purification
Dilution Range > 95% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 7040
Gene Symbol TGFB1
Uniprot ID TGFB1_HUMAN
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081