

SARS-COV-2 Nucleocapsid (B.1.1.529/Omicron) protein (Recombinant) (C-His) (STJP001478)

STJP001478

GENERAL INFORMATION





Product Type Proteins
Short Description
Host/Source E. coli

PRODUCT PROPERTIES

Concentration
Formulation Supplied as a 0.22 µg/ml filtered solution in 20mM Tris, 500mM NaCl, 0.1mM EDTA, 10% Glycerol, pH 8.0.
Purification
Dilution Range > 95% by SDS-PAGE.
Storage Instruction Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081