

Human GITR Ligand/TNFSF18 Trimer protein (Recombinant) (N-His-Flag) {Biotin} (STJP001553)

STJP001553

GENERAL INFORMATION




Product Type	Proteins
Short Description	Recombinant-Human GITR Ligand/TNFSF18 Trimer-N-His-Flag protein was developed from hek293 and has a target region of N-His-Flag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	8995
Gene Symbol	TNFSF18
Uniprot ID	TNF18_HUMAN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081