

Human IL-15RA/IL-15 R alpha/CD215 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001659)

STJP001659

GENERAL INFORMATION


Product Type	Proteins
Short Description	Recombinant-Human IL-15RA/IL-15 R alpha/CD215-C-His-Avi protein was developed from hek293 and has a target region of C-His-Avi. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES


Concentration	
Formulation	Lyophilized from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.
Instruction	

TARGET INFORMATION


Gene ID	3601
Gene Symbol	IL15RA
Uniprot ID	I15RA_HUMAN
Immunogen Sequence	



Biotinylated Human IL-15RA on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Human IL-15RA is greater than 95% as determined by SEC-HPLC.



Immoblized Human IL-15 at 0.5 µg/ml (100 µg/ml) on the plate. Dose response curve for Biotinylated Human IL-15RA, His Tag with the EC50 of 18.1 ng/ml determined by ELISA.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081