

Human LILRA2/CD85h/ILT1 protein (Recombinant) (C-His-Avi) {Biotin} (STJP001667)

STJP001667

GENERAL INFORMATION



Product Type	Proteins
Short Description	Recombinant-Human LILRA2/CD85h/ILT1-C-His-Avi protein was developed from hek293 and has a target region of C-His-Avi. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilized from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	11027
Gene Symbol	LILRA2
Uniprot ID	LIRA2_HUMAN
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081