

Human IGF2R protein (Recombinant) (C-His-Avi) {Biotin} (STJP001699)

STJP001699

GENERAL INFORMATION



Product Type	Proteins
Short Description	Recombinant-Human IGF2R-C-His-Avi protein was developed from hek293 and has a target region of C-His-Avi. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

TARGET INFORMATION

Gene ID	3482
Gene Symbol	IGF2R
Uniprot ID	MPRI_HUMAN
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081