

## Human FAP protein (Recombinant) (N-His-Avi) {Biotin} (STJP001716)

STJP001716

### GENERAL INFORMATION




<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human FAP-N-His-Avi protein was developed from hek293 and has a target region of N-His-Avi. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage Instruction</b>	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.

### TARGET INFORMATION

<b>Gene ID</b>	2191
<b>Gene Symbol</b>	FAP
<b>Uniprot ID</b>	SEPR_HUMAN
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081