

Mouse MSLN/Mesothelin protein (Recombinant) (C-His) {Biotin} (STJP001749)

STJP001749

GENERAL INFORMATION


Product Type	Proteins
Short Description	Recombinant-Mouse MSLN/Mesothelin-C-His protein was developed from hek293 and has a target region of C-His. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES


Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least 1 year from date of receipt at -20°C.
Instruction	

TARGET INFORMATION

Gene ID	56047
Gene Symbol	Msln
Uniprot ID	MSLN_MOUSE
Immunogen	
Sequence	



Biotinylated Mouse MSLN on Tris-Bis PAGE under reduced condition. The purity is greater than 95%.



The purity of Biotinylated Mouse MSLN is greater than 95% as determined by SEC-HPLC.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081