

Mouse LIMP-II protein (Recombinant) (STJP001771)

STJP001771

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse LIMP-II-protein was developed from. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration
Formulation
Purification
Dilution Range
Storage Instruction

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081