

## Human CD117 (SCFR) protein (Recombinant) (STJP001773)

STJP001773

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human CD117 (SCFR)-protein was developed from. For use in research applications.
<b>Host/Source</b>	

### PRODUCT PROPERTIES

<b>Concentration</b>
<b>Formulation</b>
<b>Purification</b>
<b>Dilution Range</b>
<b>Storage Instruction</b>

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen Sequence</b>