

HLA-A*011:01 KRAS G12V peptide complex peptide (Recombinant) {Biotin} (STJP001781)

STJP001781

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Recombinant-HLA-A*011:01 KRAS G12V peptide complex is synthetically produced from the Human HLA-A*011:01 AA25-305/B2M protein AA2-119 sequence and is suitable for use in western blot applications.
Host/Source	HEK293 cells NA

PRODUCT PROPERTIES

Concentration	
Formulation	Lyophilised from 0.2 µm filtered PBS solution pH7.4 NA
Purification	
Dilution Range	purity was determined to be >95% by nonreducing PAGE and HPLC (see below methods) NA
Storage Instruction	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C. NA

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence