

## HLA-A\*02:01 NY-ESO-1 peptide complex peptide (Recombinant) (STJP001785)

STJP001785

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Recombinant-HLA-A*02:01 NY-ESO-1 peptide complex is synthetically produced from the Human HLA-A*02:01 AA25-305/B2M protein AA2-119 sequence and is suitable for use in western blot applications.
<b>Host/Source</b>	HEK293 cells NA

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Lyophilised from 0.2 µm filtered PBS solution pH7.4 NA
<b>Purification</b>	
<b>Dilution Range</b>	purity was determined to be >95% by nonreducing PAGE and HPLC (see below methods) NA
<b>Storage Instruction</b>	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C C. NA

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081