

HLA-E*01:03 peptide complex peptide (Recombinant) (STJP001796)

STJP001796

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Recombinant-HLA-E*01:03 peptide complex is synthetically produced from the sequence and is suitable for use in western blot applications.
Host/Source	

PRODUCT PROPERTIES

Concentration
Formulation
Purification
Dilution Range
Storage Instruction

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081