

## HLA-E\*01:03 peptide complex Tetramer peptide (Recombinant) (STJP001797)

STJP001797

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Recombinant-HLA-E*01:03 peptide complex Tetramer is synthetically produced from the sequence and is suitable for use in western blot applications.
<b>Host/Source</b>	

### PRODUCT PROPERTIES

Concentration  
Formulation  
Purification  
Dilution Range  
Storage  
Instruction

### TARGET INFORMATION

Gene ID  
Gene Symbol  
Uniprot ID  
Immunogen  
Sequence